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Instructions for Use EN

GOeasy Basic Bond

One component light cure adhesive for total-etch technique

GOeasy Basic Bond is a simple to use light cure one component bonding agent. It is designed for strong bonding of light cure composites and compomers to etched enamel and dentine and to nonprecious and precious metals.

Its strong adhesion to etched enamel or dentine proceeds on principles similar to that occurring with glassionomer cements. Good, long lasting adhesive strength and good biocompatibility are attained by polycarboxylic acids.

GOeasy Basic Bond is compatible with all current brands of visible light cure composite materials. It is ethanol-based and hydrophilic. GOeasy Basic Bond can be used on slightly moist dentine surfaces (wet bonding technique).

Contraindications

Do not use GOeasy Basic Bond if the recommended working technique is not possible or if the patient is known to be allergic to any of the ingredients.

Side effects

In singular cases, GOeasy Basic Bond may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin

layer of calcium hydroxide material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of GOeasy Basic Bond. Neither store the material in proximity of eugenol containing products, nor let the material allow coming into contact with materials containing eugenol.

Application

1. Isolation

Rubber dam is the recommended method of isolation.

2. Cavity preparation

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (**0.5 - 1.0 mm**) bevel placed in the enamel to increase the surface area for greater bond strength.

3. Pulp protection

Cavity floor of deep excavations should be covered with a thin layer of calcium hydroxide material (e.g. **GOeasy Ca-Line**).

4. Enamel and dentine conditioning

Recommended is the total etch technique. Apply **GOeasy Etch** onto the enamel and dentine surfaces beginning with the enamel bevels. Leave the **GOeasy Etch** in place for **20 seconds**. On primary teeth and teeth high in fluoride, enamel etching up to 60 seconds is possible. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream, but do not desiccate. A slightly wet dentine surface

is important for the function of GOeasy Basic Bond. The etched enamel bevel should have a chalky white appearance.

Etching precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

For the typical indications of compomers etching is mostly not recommended, but it improves considerably the bonding with tooth material.

5. Application of GOeasy Basic Bond

Apply GOeasy Basic Bond generously with a brush onto the enamel and dentine surfaces for **30 seconds** with agitation. The material should build a homogeneous layer. Remove excess material carefully. Dry cautiously with oil free air for about 15 seconds to remove all volatile components and to disperse the adhesive to an even layer. Do not desiccate the dentine.

Cure the GOeasy Basic Bond coating by exposing its entire area to a dental halogen light unit for **20 seconds** before application of a second layer of GOeasy Basic Bond.

Notes: Do not rinse off the GOeasy Basic Bond! If not used immediately, place dispensed GOeasy Basic Bond in subdued light to prevent premature polymerization by

incident light. The GOeasy Basic Bond will not selfcure.

6. Application of a second layer of GOeasy Basic Bond

Apply again the GOeasy Basic Bond generously with a brush onto the adhesive surfaces as described under **point 5**. before placement of a light cure composite.

7. Restorative placement

Apply the restorative material according to the instructions of the manufacturer. Best results are obtained with application of a thin layer of a flowable composite (e.g. **GOeasy Flow**) followed by the application of a moldable composite (e.g. **GOeasy Fill**). Light cure each composite layer separately according to the corresponding user instructions.

Warnings

- Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

Storage

Do not store above 25 °C (77 °F)! Do not use after expiry date. Avoid storage in direct sunlight.

Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental professional use only!

The instructions for use has to be kept for the duration of the application.

Caution:

Federal law restricts this device to sale by or on the order of a dentist.