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## Instructions for Use EN

### GOeasy Etchbond Duo

#### Two component self-etching dual cure adhesive system

**GOeasy Etchbond Duo** is a simple to use self-etching dual cure adhesive. It is designed for strong bonding of composites, compomers and resin modified glass ionomer cements to enamel and dentine. GOeasy Etchbond Duo works optimally under slightly wet conditions, which are typical during clinical procedures. GOeasy Etchbond Duo consists of the components PART A and PART B that were mixed before application. It is bonding to dual cure, self-cure and light cure composite restorative materials.

#### Contraindications

Do not use GOeasy Etchbond Duo if the recommended working technique is not possible or if the patient is known to be allergic to any of the ingredients.

#### Side effects

In singular cases, GOeasy Etchbond Duo may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of calcium hydroxide material.

#### Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of GOeasy Etchbond Duo. Neither store the material in proximity of eugenol containing products, nor let the material allow coming into contact with materials containing eugenol.

#### Application

##### 1. Isolation

Rubber dam is the recommended method of isolation.

##### 2. Cavity preparation

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (**0.5 - 1.0 mm**) bevel placed in the enamel to increase the surface area for greater bond strength.

##### 3. Pulp protection

Cavity floor of deep excavations should be covered with a thin layer of calcium hydroxide material (e.g. GOeasy Ca-Line).

##### 4. Application of GOeasy Etchbond Duo

One drop of GOeasy Etchbond Duo PART A and one drop of GOeasy Etchbond Duo PART B were combined in a mixing pallet and mixed for 5-10 seconds.

Apply the homogeneous mixture generously with a brush onto the slightly wet enamel and dentine surfaces (as they are usual in clinical

procedures) for 30 seconds with agitation. The material should build a homogeneous layer. Air thin for 10 seconds to remove the volatile components and to disperse the adhesive. Then light cure for 20 sec with a dental halogen light unit or an LED (wavelength 400–500 nm, light intensity min. 1000 mW/cm<sup>2</sup>) and place the restorative material.

**For a maximal adhesion it is strongly recommended to use the light cure mode.**

If light cure is absolutely impossible, the adhesive will also cure in the autocure mode. After application of the adhesive air thin gently to remove all volatile components. There must remain a sticky layer. Then apply the restorative material directly.

#### Notes:

Do not interchange lids of the bottles, because this can lead to a cross-contamination of the liquids.

##### 5. Restorative placement

Apply the restorative material according to the instructions of the manufacturer.

Best results are obtained with application of a thin layer of a flowable composite (e.g. **GOeasy Flow**) followed by the application of a moldable composite (e.g. **GOeasy Fill**). Light cure each composite layer separately according to the corresponding user instructions.

#### Warnings

- Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

#### Storage

Do not store above 25 °C (77 °F)!  
Avoid storage in direct sunlight.  
Do not use after expiry date.

#### Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture.  
P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

#### Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

#### Keep away from children!

#### For dental professional use only!

The instructions for use has to be kept for the duration of the application.

#### Caution:

**Federal law restricts this device to sale by or on the order of a dentist.**