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### Instructions for Use EN

## **GOeasy Glass Cem Classic** **Conventional glass ionomer luting** **cement**

**GOeasy Glass Cem Classic** is a fluoride containing radiopaque classical strontium glass ionomer cement designed for cementing a wide variety of prosthetic appliances and materials. As a glass ionomer cement, it features good chemical adhesion and biocompatibility. It bonds to enamel, dentine, non-precious metals and porcelain.

GOeasy Glass Cem Classic meets the requirements of ISO 9917-1.

### **Indications**

- cementing of crowns, bridges, inlays and onlays (all types: metal, resin, ceramic fused to metal and ceramic)
- cementing of endodontic posts and orthodontic bands
- liner under composite fillings

### **1. Tooth Preparation**

Prepare tooth in usual manner. Clean prepared tooth with pumice and water only. Rinse thoroughly and dry, but do not desiccate.

By using a cotton pellet apply a dentine conditioner according to the corresponding instructions to remove the smear layer. Rinse the dentine conditioner with water and dry it in an airstream, but do not desiccate.

Pulp capping with GOeasy Glass Cem Classic is contra indicated.

To deep areas of possible pulpal exposure apply a calcium hydroxide liner (e.g. GOeasy Ca-Line).

### **2. Dispensing and Mixing**

The powder/liquid ratio to achieve a suitable consistency is 1.8 / 1.0. This can be obtained by mixing of one scoop

(blue colored) of powder and two drops of liquid.

For accurate dispensing of GOeasy Cem Classic Powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder by using the collar at the top of the bottle and carry it onto the mixing pad. Avoid compressing powder into the spoon with the inside wall of the bottle.

For dispensing of GOeasy Cem Classic Liquid turn the bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense the drops. If any bubbles are present, lightly tap the bottle with the fingers holding it. Discount drops that contain bubbles and are obviously not full-sized.

Use a small spatula to rapidly mix half the powder into the liquid. Then mix the remaining powder into the mixed paste to achieve a glossy appearance. Total mixing time is **30 seconds**.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

### **3. Cementation Technique**

Mix the required amount of GOeasy Glass Cem Classic. Coat the prepared bonding surface of the restoration with **1 mm** of cement (**Do not overfill!**) and seat immediately (**working time 1:25 min** from start of mixing at 23°C). Seat the restoration by using steady, even pressure. Net setting time is about **4:30 min**.

Remove excess cement at the first setting stage.

Maintain isolation until the set of the cement is verified (approx. 4:30 min).

**Note:** Higher temperatures will shorten the working time, lower temperatures will prolong the working time. An overextended working time will cause the loss of adhesion to the dental enamel and the dentin.

### **4. Conclusive Notes**

This product is to be applied only by a dental professional in the way as described in the instructions.

Do not use GOeasy Glass Cem Classic with patients who show an allergy to the material. In case of allergic reactions immediately stop the application, and advise the patient to consult a physician. An operator, who has a history of allergy to glass ionomer-cements should not handle GOeasy Glass Cem Classic.

Do not allow the liquid or cement mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.

Avoid eye contact of the liquid or cement mixture. In case of contact, immediately flush with water and seek medical treatment.

Do not mix the powder or liquid of GOeasy Glass Cem Classic with any other glass-ionomer product.

### **5. Storage**

Store GOeasy Glass Cem Classic in a cool and dark place at 4-25 °C (39-77 °F). Temperature should not exceed 25 °C (77 °F). Do not use after expiry date.

### **Warranty**

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L.

Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

### **Limitation of Liability**

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

### **Keep away from children!**

### **For dental professional use only!**

The instructions for use has to be kept for the duration of the application.

### **Caution:**

**Federal law restricts this device to sale by or on the order of a dentist.**