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### GOeasy Glass Core & Line

#### Core build-up glass ionomer cement

**GOeasy Glass Core & Line** is a radiopaque fluoride containing glass ionomer cement. The excellent chemical bonding to tooth material (dentine and enamel) requires no enamel and dentine etching and serves for excellent marginal seal. The early resistance to water uptake permits a one-visit-treatment. It is strongly radiopaque, thus is ensuring easy postoperative diagnosis. Due to its biocompatibility, pulp irritations are minimized.

GOeasy Glass Core & Line meets the requirements of ISO 9917-1.

#### Indications

- Core build-ups
- Intermediate restorative and sandwich material for heavy stress bearing Class I and Class II cavities

#### 1. Tooth Preparation

Clean the cavity preparation with pumice and water. Rinse thoroughly and dry, but do not desiccate.

By using a cotton pellet apply a dentine conditioner according to the corresponding instructions to remove the smear layer. Rinse the dentine conditioner with water and dry it in an airstream, but do not desiccate.

Pulp capping with GOeasy Glass Core & Line is contraindicated. To deep areas of possible pulpal exposure apply a calcium hydroxide liner (e.g. **GOeasy Ca-Line**).

#### 2. Dispensing and Mixing

The powder/liquid ratio to achieve a suitable consistency is 3.6 / 1.0. This can be obtained by mixing of 1 level scoop (blue) of powder and 1 drop of liquid.

A condensable thick mix is required to minimize abrasive loss, but a glossy surface is also necessary for strong chemical bond strength. For application the appropriate powder / liquid ratio is strongly recommended.

Shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad provided. Avoid compressing powder into the spoon with the inside wall of the bottle.

Turn the liquid bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. Discard drops that are obviously not full-sized.

Use a small spatula to rapidly mix half the powder into the liquid. Then mix the remaining powder into the mixed paste to achieve a glossy appearance. Total mixing time is **30 sec**.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

#### 3. Filling

Remove moisture from the cavity surface either with a cotton pellet, or very gently with an air syringe blast. It is essential not to desiccate the dentine and enamel surfaces.

Mix the required amount of GOeasy Glass Core & Line and apply it by using a suitable placement instrument. Please see to it, that no air bubbles will be incorporated.

Form the material with a placement or forming instrument and use a matrix.

**Total working time from beginning the mixing is 1:30 min (at 23°C).**

#### Note:

Higher temperatures will shorten the working time; lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the dental enamel and the dentine.

Remove the matrix approx. **3:30 minutes** after having finished the application (**net setting time**). **Total setting time from start of mixing is approx. 5:00 min.**

#### 4. Finishing

Final finishing and polishing can begin from about 6 minutes after start of mixing. Smooth with fine diamond burs, then polish with finishing and polishing discs with graded grain size.

Apply a final coating of varnish to the final finished surface of the restoration.

Instruct the patient not to expose the restoration to any pressure for one hour.

#### 5. Conclusive Notes

This product is to be applied only by a dental professional in the way as described in this instruction.

Do not use GOeasy Glass Core & Line with patients who show an allergy to the material. In case of allergic reactions immediately stop the application, and advise the patient to consult a physician. An operator, who has a history of allergy to glass ionomer cements should not handle GOeasy Glass Core & Line.

Do not allow the liquid or cement mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.

Avoid eye contact of the liquid or cement mixture.

In case of contact, immediately flush with water and seek medical treatment.

Do not mix the powder or liquid of GOeasy Glass Core & Line with any other glass ionomer product.

#### 6. Storage

Store GOeasy Glass Core & Line in a cool place at 4-25 °C (39-77 °F). Temperature should not exceed 25 °C (77 °F). Do not use after expiry date.

#### Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product

for user's application. If this product is defective within the warranty period, your exclusive remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

#### Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

#### Keep away from children!

#### For dental professional use only!

The instructions for use has to be kept for the duration of the application.

#### Caution:

**Federal law restricts this device to sale by or on the order of a dentist.**