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Instructions for Use EN

GOeasy Glass Fill Extra Esthetic conventional glass ionomer filling cement

GOeasy Glass Fill Extra is a classical restorative glass ionomer cement with improved esthetic and good physical properties. Besides its excellent fluoride release and excellent biocompatibility as glass ionomer cement GOeasy Glass Fill Extra also has good chemical bonding to dentin and enamel and a tight seal at the dentinal margins. Therefore it requires no enamel and dentin etching. Because of its radiopacity it ensures easy postoperative diagnosis.

GOeasy Glass Fill Extra meets the requirements of ISO 9917-1.

Indications

- linings for composite fillings
- primary tooth fillings
- core build-ups
- dental neck fillings
- uniplanar and multiplanar temporary fillings
- small class I fillings

1. Tooth Preparation

Clean the cavity preparation with pumice and water. Rinse thoroughly and dry, but do not desiccate.

If desired, place a matrix band.

By using a cotton pellet apply a dentine conditioner according to the corresponding instructions to remove the smear layer. Rinse the dentine conditioner with water and dry it in an airstream, but do not desiccate.

Pulp capping with GOeasy Glass Fill Extra is contra indicated. To deep areas of possible pulpal exposure apply a calcium hydroxide liner (e.g. GOeasy Ca-Line).

2. Dispensing and Mixing

The powder/liquid ratio to achieve a suitable consistency is 3.6 / 1.0. This can be obtained by mixing of one level scoop (blue colored) of powder and one drop of liquid.

A condensable thick mix is required to minimize abrasive loss, but a glossy surface is also need for strong chemical bond strength. Using the appropriate powder / liquid ratio is therefore strongly recommended.

Shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad provided. Avoid compressing powder into the spoon with the inside wall of the bottle.

Turn the liquid bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. Discard drops that are obviously not full-sized.

Use a small spatula to rapidly mix all the cement powder into the liquid. Usually one scoop powder/ one drop liquid should provide for sufficient amount of mixed cement. The mixed cement should be thixotropic and have a smooth consistency and glossy appearance. Total mixing time is **30 seconds**.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

3. Filling

Apply the mixed material into the prepared cavity by using a suitable placement instrument.

Please see to it, that no air bubbles will be incorporated.

Form the cement with a placement or forming instrument within the **working time (1:15 minutes at 23 °C or 74°F)** and use a matrix. Net setting time is about **3:30 minutes**.

Immediately after setting, a varnish should be placed onto the surface of the restoration.

Note:

Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the dental enamel and the dentin.

Remove the matrix when the cement has achieved clinical set (**3:30 minutes** after application)

4. Finishing

Final finishing and polishing can begin from about 6 minutes after start of mixing. Smooth with fine diamond burs, then polish with finishing and polishing discs with graded grain size.

Spray preparation dust away with water and dry the surface in an oilfree airstream.

Apply a final coating of varnish to the final finished surface of the restoration and instruct the patient not to expose it to any pressure for one hour.

5. Conclusive Notes

This product is to be applied only by a dental professional in the way as described in the instructions.

Do not use GOeasy Glass Fill Extra with patients who show an allergy to the material. In case of allergic reactions

immediately stop the application, and advise the patient to consult a physician. An operator, who has a history of allergy to glass ionomer cements should not handle GOeasy Glass Fill Extra.

Do not allow the liquid or cement mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.

Avoid eye contact of the liquid or cement mixture. In case of contact, immediately flush with water and seek medical treatment.

Do not mix the powder or liquid of GOeasy Glass Fill Extra with any other glass ionomer product.

6. Storage

Store GOeasy Glass Fill Extra in a cool place at 4-25 °C (39-77 °F). Temperature should not exceed 25 °C (77 °F). Do not use after expiry date.

Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect,

special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental professional use only!

The instructions for use has to be kept for the duration of the application.

Caution:

Federal law restricts this device to sale by or on the order of a dentist.