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Light cure resin modified glass ionomer filling cement

GOeasy ResinGlass Fill is a light cure resin reinforced glass ionomer filling cement. Beside its excellent fluoride release and the excellent biocompatibility as glass ionomer cement GOeasy ResinGlass Fill has also good chemical bonding to dentine and enamel and a tight seal at the dentinal margins. Because of its radiopacity it ensures easy postoperative diagnosis.

GOeasy ResinGlass Fill meets the requirements of ISO 9917-2.

Indications

- linings for composite fillings
- primary tooth fillings
- core build-ups
- dental neck fillings
- uniplanar and multiplanar temporary fillings
- small class I and II fillings

1. Tooth Preparation

Clean the cavity preparation with pumice and water. Rinse thoroughly and dry gently, but do not desiccate.

Apply **GOeasy ResinGlass Prime** with a suitable brush onto the moist enamel and dentine surfaces for 30 seconds with agitation. The material should build a homogeneous layer. Air thin gently for 10 seconds to remove the volatile components and to disperse the adhesive. Then light cure with a suitable dental light unit for 10 seconds.

Pulp capping with GOeasy ResinGlass Fill is contraindicated. To deep areas of

possible pulpal exposure apply a small amount of a calcium hydroxide liner (e.g. **GOeasy Ca-Line**)

2. Dispensing and Mixing

The powder/liquid ratio is 2.6/1.0. This can be obtained by mixing 1 level (blue) scoop of powder and 2 drops of liquid.

For accurate dispensing of GOeasy ResinGlass Fill powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad. Avoid compressing powder into the spoon with the inside wall of the bottle.

For dispensing of GOeasy ResinGlass Fill liquid turn the bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. Discard drops that contain bubbles and are obviously not full-sized.

Use a small spatula to rapidly mix all the cement powder into the liquid. The mixed cement should be thixotropic and have a smooth consistency and glossy appearance. Total mixing time is **30 seconds**.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

3. Filling and Finishing

Mix the required amount of GOeasy ResinGlass Fill and apply it into the prepared cavity by using a suitable

placement instrument.

Please see to it, that no air bubbles will be incorporated.

GOeasy ResinGlass Fill should be placed in the cavity within the working time (**2:00 minutes** from start of mixing at 23 °C or 74°F). A transparent matrix may be used.

Note: Higher temperatures will shorten the working time, lower temperatures will prolong the working time. An overextended working time will cause the loss of adhesion to the dental enamel and the dentine.

Cure **for 20 sec** with a visible light-curing device. Net setting time without any light is about **4:30 min**.

If the cavity is more than **2 mm** deep, curing in thin incremental layers is strongly recommended.

After GOeasy ResinGlass Fill has set, remove matrix and perform finishing under water spray using standard techniques.

4. Conclusive Notes

This product is to be applied only by a dental professional in the manner as described in this instruction.

Do not use GOeasy ResinGlass Fill with patients who show an allergy to the material. In case of allergic reactions immediately stop the application, and advise the patient to consult a physician. An operator, who has a history of allergy to glass ionomer cements, should not

handle GOeasy ResinGlass Fill.

Do not allow the cement mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.

Avoid eye contact of the cement mixture. In case of contact, immediately flush with water and seek medical treatment.

Do not mix the powder or liquid of GOeasy ResinGlass Fill with any other glass ionomer product.

Direct contact with eugenol-based products is to be avoided. Eugenol inhibits the setting of LC glass ionomer cements.

5. Storage

Store GOeasy ResinGlass Fill in a cool and dark place at 4-25 °C (39-77 °F). Temperature should not exceed 25 °C (77 °F). Do not use after expiry date.

Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental professional use only!

The instructions for use has to be kept for the duration of the application.

Caution:

Federal law restricts this device to sale by or on the order of a dentist.